


Product Datasheet

ECFP-Tag antibody (orb763636)

Description	Mouse polyclonal antibody to ECFP-Tag
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Immunogen	Recombinant Protein
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Source	Mouse
Dilution Range	WB 1:2,000-5,000 IF 1:200 IHC 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of ECFP Recombinan...