

Product Datasheet

EYFP-Tag antibody (orb763647)

Description	Rabbit polyclonal antibody to EYFP-Tag
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB 1:3,000-5,000 IHC 1:200 IF 1:200 IP 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of EYFP recombinan...