

Product Datasheet

IL17/IL17A Antibody (orb76367)

| | |
|----------------------------|---|
| Description | IL17/IL17A Antibody |
| Species/Host | Rabbit |
| Reactivity | Mouse |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | E. coli-derived mouse IL-17 recombinant protein (Position: A26-A158). |
| Form/Appearance | Lyophilized |
| Concentration | Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml. |
| Storage | Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles. |
| Note | For research use only |
| Application notes | Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml. |
| Isotype | Rabbit IgG |
| Clonality | Polyclonal |
| MW | 17490 MW |
| Uniprot ID | Q62386 |
| Expiration Date | 12 months from date of receipt. |



Figure.
Western blot
analysis of
IL-17 u...