

## Product Datasheet

### Histone H3 (phospho-Thr32) antibody (orb763674)

**Description** Rabbit polyclonal antibody to Histone H3 (phospho-Thr32)

**Species/Host** Rabbit

**Reactivity** Human, Mouse, Rat

**Conjugation** Unconjugated

**Tested Applications** IHC, WB

**Immunogen** Synthetic Peptide

**Preservatives** PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

**Concentration** 1mg/ml

**Storage** Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note** For research use only

**Isotype** IgG


**Clonality** Polyclonal

**MW** 15-17kDa


**Uniprot ID** **P68431**

**Dilution Range** WB 1:1,000-2,000 IHC 1:200-500

**Expiration Date** 12 months from date of receipt.



Western blot analysis of extracts from H...



Immunohistochemical analysis of paraffin...