



Product Datasheet

I kappa B alpha antibody (orb763735)

Description	Rabbit polyclonal antibody to I kappa B alpha
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant Protein
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40kDa
Uniprot ID	P25963
Dilution Range	WB 1:2,000-5,000 IHC 1:200-500
Expiration Date	12 months from date of receipt.



Western blot analysis of 1) HeLa, 2) PC1...



Immunohistochemical analysis of paraffin...