


Product Datasheet

SARS-E2 Antibody (orb76374)

Description	Mouse monoclonal antibody to SARS-E2
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified recombinant fragment of SARS-E2 glycoprotein precursor expressed in E. Coli.
Preservatives	Subclonal supernatant.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Dilution Range	ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of SARS-E2GP3 recombinant protein.