


## Product Datasheet

### SARS-M Antibody (orb76375)

|                            |                                                                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to SARS-M                                                                                                     |
| <b>Species/Host</b>        | Mouse                                                                                                                                   |
| <b>Reactivity</b>          | Virus                                                                                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                                                                            |
| <b>Tested Applications</b> | ELISA                                                                                                                                   |
| <b>Immunogen</b>           | Purified recombinant fragment of SARS-m protein expressed in E. Coli.                                                                   |
| <b>Preservatives</b>       | Subclonal supernatant.                                                                                                                  |
| <b>Concentration</b>       | 1mg/ml                                                                                                                                  |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                   |
| <b>Isotype</b>             | Mouse IgG1                                                                                                                              |
| <b>Clonality</b>           | Monoclonal                                                                                                                              |
| <b>Dilution Range</b>      | ELISA: 1/10000                                                                                                                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                         |



Western blot analysis of SARS-m protein expressed in E. Coli.