



## Product Datasheet

### Cav1.3 antibody (orb763785)

<b>Description</b>	Rabbit polyclonal antibody to Cav1.3
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC
<b>Immunogen</b>	Synthetic Peptide
<b>Preservatives</b>	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Store at -20°C. Avoid repeated freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	245kDa
<b>Uniprot ID</b>	<a href="https://www.uniprot.org/entry/Q01668">Q01668</a>
<b>Dilution Range</b>	IHC 1:100-200
<b>Expiration Date</b>	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...