


---

## Product Datasheet

### Cleaved-Caspase-4 p20 (Q81) antibody (orb763894)

**Description**


Rabbit polyclonal antibody to Cleaved-Caspase-4

Western Blot analysis of various cells u...



Western Blot analysis of 3T3 cells using...



Western blot analysis of lysates from Ju...

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse

**Conjugation**

Unconjugated

**Tested Applications**

ELISA, IHC-P, WB

**Immunogen**

The antiserum was produced against synthesized peptide derived from human Caspase 4. AA range:62-111

**Preservatives**

PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

**Concentration**

1 mg/ml

**Storage**

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**MW**

43+22

**Expiration Date**

12 months from date of receipt.