

Product Datasheet

Myc-Tag antibody (orb764018)

Description	Rabbit polyclonal antibody to Myc-Tag
Species/Host	Rabbit
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide EQKLISEEDL coupled to KLH.
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.