

Product Datasheet

Alliinase antibody (orb764035)

| | |
|----------------------------|--|
| Description | Rabbit polyclonal antibody to Alliinase |
| Species/Host | Rabbit |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | Purified Protein |
| Preservatives | PBS with 0.02% sodium azide and 50% glycerol pH 7.4. |
| Storage | Store at -20°C. Avoid repeated freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 50 |
| Expiration Date | 12 months from date of receipt. |