

Product Datasheet

Histone H3 (phospho-Thr45) antibody (orb764069)

Description	Rabbit polyclonal antibody to Histone H3 (phospho-Thr45)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide of Histone H3 (Phospho Thr45)
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	15-17
Expiration Date	12 months from date of receipt.