



## Product Datasheet

# NF kappa B-p65 (phospho-Ser529) antibody (orb764245)


<b>Description</b>	Rabbit polyclonal antibody to NF kappa B-p65 (phospho-Ser529)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Monkey, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IF, IHC-P, WB
<b>Immunogen</b>	Synthesized phospho-peptide around the phosphorylation site of human NFκB-p65 (phospho Ser529)
<b>Preservatives</b>	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	60
<b>Expiration Date</b>	12 months from date of receipt.



Immunofluorescence analysis of human-lun...



Immunofluorescence analysis of human-lun...



Immunofluorescence analysis of rat-heart...