



Product Datasheet

Sp1 (phospho-Thr739) antibody (orb764278)


Description	Rabbit polyclonal antibody to Sp1 (phospho-Thr739)
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human SP1 around the phosphorylation site of Thr739. AA range:706-755
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	81
Expiration Date	12 months from date of receipt.



Western Blot analysis of various cells
u...



Enzyme-Linked Immunosorbent Assay (Phosp...



Western blot analysis of lysates from CO...