


Product Datasheet


NOS3 (phospho-Ser1177) antibody (orb764317)

Description	Rabbit polyclonal antibody to NOS3 (phospho-Ser1177)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human eNOS around the phosphorylation site of Ser1176. AA range:1144-1193
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	133
Expiration Date	12 months from date of receipt.




Immunofluorescence analysis of A549. 1,

...



Immunofluorescence analysis of HeLa cell...



Western Blot analysis of various cells u...