


Product Datasheet


14-3-3 zeta (phospho-Ser58) antibody (orb764378)

Description	Rabbit polyclonal antibody to 14-3-3 zeta (phospho-Ser58)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human 14-3-3 zeta around the phosphorylation site of Ser58. AA range:24-73
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	30
Expiration Date	12 months from date of receipt.




Western Blot analysis of various cells

U...



Western Blot analysis of Jurkat cells us...



Enzyme-Linked Immunosorbent Assay (Phosp...