

## Product Datasheet

# Amyloid-beta (phospho-Thr743) antibody (orb764380)

**Description** Rabbit polyclonal antibody to Amyloid-beta (phospho-Thr743)

**Species/Host** Rabbit

**Reactivity** Human, Mouse, Rat

**Conjugation** Unconjugated

**Tested Applications** ELISA, IF, IHC-P, WB

**Immunogen** The antiserum was produced against synthesized peptide derived from human Amyloid beta A4 around the phosphorylation site of Thr743/668. AA range:711-760

**Preservatives** PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

**Concentration** 1 mg/ml

**Storage** Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.


**Note** For research use only

**Isotype** IgG


**Clonality** Polyclonal

**MW** 140


**Expiration Date** 12 months from date of receipt.



Western Blot analysis of various cells u...



Immunohistochemical analysis of paraffin...



Enzyme-Linked Immunosorbent Assay (Phosp...