



Product Datasheet

Cytokeratin 17 antibody (orb765002)


Description	Rabbit polyclonal antibody to Cytokeratin 17
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human Keratin 17. AA range:381-430
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	48
Expiration Date	12 months from date of receipt.



Western Blot analysis
of various cells u...



Immunofluorescence
analysis of HepG2
cel...



Immunohistochemistry
analysis of paraffi...