



Product Datasheet

JAB1 antibody (orb765529)


Description	Rabbit polyclonal antibody to JAB1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human COPS5. AA range:161-210
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	38
Expiration Date	12 months from date of receipt.



Western Blot analysis of K562 cells using JAB1 antibody.



Western Blot analysis of various cells using JAB1 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using JAB1 antibody.