



Product Datasheet

PKM2(phospho-Ser37) Antibody (orb76553)

Description	PKM2(phospho-Ser37) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P14618
Expiration Date	12 months from date of receipt.



Western blot analysis of PKM2(phospho-Se...



Immunofluorescence staining of methanol-...