


Product Datasheet

Interferon- α Receptor Type I Subunit I (phospho-Ser535/Ser539) Antibody (orb76554)

Description	Interferon- α Receptor Type I Subunit I (phospho-Ser535/Ser539) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P17181
Expiration Date	12 months from date of receipt.



Western blot analysis of Interferon- α Re...