


Product Datasheet

PTH rP antibody (orb76579)

Description	Mouse monoclonal antibody to PTH rP
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC-P, WB
Immunogen	23 amino acid peptide (corresponding to amino acids 146-169) of the hPTHrP receptor
Target	Parathyroid Hormone Receptor 1
Preservatives	Purified antibody containing PBS 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	P3/X63.Ag8.653 mouse myeloma cells with spleen cells from a BALB/c;
Isotype	Ig
Clonality	Monoclonal
Uniprot ID	p12272
Entrez	5744
Dilution Range	Optimal antibody dilution should be determined by titration
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of PTHLH an...



Line graph illustrates about the Ag-Ab r...




Figure 1. Western blot using anti-3D1.1 immunant (SN) or purified antibody (Pur) with extracts from HEK cells or HEK cells transfected with hPTH1r/PTHrP receptor cDNA (C21).

Western blot analysis of HEK cells using...