


Product Datasheet

CCR9 antibody (orb76597)

Description	Rat monoclonal antibody to CCR9
Species/Host	Rat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-Fr
Immunogen	Synthetic peptide: PTDFTSIPNMADDYGSEST, corresponding to N terminal amino acids 3-22 of Mouse CCR9
Target	CCR9
Preservatives	Tissue culture supernatant in RPMI 1640 5 % FCS plus 0.1% sodium azide;
Concentration	20ug/ml
Note	For research use only
Application notes	P3-x63Ag.8.653;
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	Q63190
Entrez	10803 , 12769
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try a dilution of 1:20
Expiration Date	12 months from date of receipt.



Confocal immunofluorescence analysis of ...