



Product Datasheet

Ribosomal Protein L17 antibody (orb766233)


Description	Rabbit polyclonal antibody to Ribosomal Protein L17
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human RPL17. AA range:101-150
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	24
Expiration Date	12 months from date of receipt.



Western Blot analysis of various cells u...



Western Blot analysis of customer's (cat...



Immunohistochemistry analysis of paraffi...