



Product Datasheet

Snk antibody (orb766347)


| | |
|----------------------------|--|
| Description | Rabbit polyclonal antibody to Snk |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | The antiserum was produced against synthesized peptide derived from human PLK2. AA range:291-340 |
| Preservatives | PBS with 0.02% sodium azide and 50% glycerol pH 7.4. |
| Concentration | 1 mg/ml |
| Storage | Store at -20°C. Avoid repeated freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 78 |
| Expiration Date | 12 months from date of receipt. |



Western Blot analysis of various cells u...



Western Blot analysis of MCF7 cells usin...



Western blot analysis of lysates from He...