


Product Datasheet

Glutathione S-transferase antibody (orb76636)

Description	Mouse monoclonal antibody to Glutathione S-transferase
Species/Host	Mouse
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	GST-fusion protein
Target	GST
Preservatives	Purified antibody (from supernatant) containing PBS + 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	X63-Ag8.653
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P08515
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline start at 1:2,000 for western blot
Expiration Date	12 months from date of receipt.



Western blot analysis of GST fusion prot...