


Product Datasheet

HIV1 Vif antibody (orb76655)

Description	Mouse monoclonal antibody to HIV1 Vif
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Bacterially expressed, hexahistidine amino-terminal tagged HIV-1 Vif protein (clade B, HXB-3 isolate)
Target	HIV1 Vif
Preservatives	Affinity purified antiserum containing PBS 0.1% sodium azide;
Concentration	1mg/ml
Note	For research use only
Application notes	Sp2/0 Ag.14
Isotype	IgG1 kappa
Clonality	Monoclonal
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try; WB: 1/1000. Predicted molecular weight: 23 kDa
Expiration Date	12 months from date of receipt.



Western blot analysis of semi purified H...