

Product Datasheet

HIV2 Nef antibody (orb76656)

Description	Rat monoclonal antibody to HIV2 Nef
Species/Host	Rat
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	TrpE-Fusion protein
Target	HIV2 Nef
Preservatives	Purified antibody containing PBS 0.1% sodium azide;
Concentration	1mg/ml
Note	For research use only
Application notes	X63Ag8.653
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	Optimal antibody dilution should be determined by titration, however try; WB: Use at 1:1000. Predicted molecular weight: 35 kDa.
Expiration Date	12 months from date of receipt.