



Product Datasheet

Stathmin-2 antibody (orb766681)


Description	Rabbit polyclonal antibody to Stathmin-2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthesized peptide derived from the Internal region of human Stathmin-2.
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20
Expiration Date	12 months from date of receipt.



Western Blot analysis of extracts from J...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...