



Product Datasheet

Prothrombin antibody (orb766731)


| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to Prothrombin |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IHC-P, WB |
| Immunogen | Synthesized peptide derived from the Internal region of human Prothrombin. |
| Preservatives | PBS with 0.02% sodium azide and 50% glycerol pH 7.4. |
| Concentration | 1 mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 70 |
| Expiration Date | 12 months from date of receipt. |



Western Blot analysis of Jurkat cells us...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...