



Product Datasheet

GGT1 antibody (orb766871)


Description	Rabbit polyclonal antibody to GGT1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the N-terminal region of human GGT1. AA range:21-70
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	61
Expiration Date	12 months from date of receipt.



Western Blot analysis of HepG2 cells using GGT1 antibody.



Immunohistochemical analysis of paraffin-embedded tissue using GGT1 antibody.



Immunohistochemical analysis of paraffin-embedded tissue using GGT1 antibody.