



Product Datasheet

IL13R alpha 2 antibody (orb767221)


Description	Rabbit polyclonal antibody to IL13R alpha 2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthesized peptide derived from Interleukin-13 receptor subunit alpha-2 at AA range: 241-290
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	44
Expiration Date	12 months from date of receipt.



Western Blot analysis of HepG2 (UV), L92...



Western Blot analysis of HEPG2-UV cells ...



Immunohistochemical analysis of paraffin...