

Product Datasheet

Cdk7 (phospho-Thr170) antibody (orb767317)

Description Rabbit polyclonal antibody to Cdk7 (phospho-Thr170)

Species/Host Rabbit

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications ELISA, IHC-P, WB

Immunogen The antiserum was produced against synthesized peptide derived from human CDK7 around the phosphorylation site of Thr170. AA range:136-185

Preservatives PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Concentration 1 mg/ml

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.


Note For research use only

Isotype IgG


Clonality Polyclonal

MW 40

Expiration Date 12 months from date of receipt.



Immunohistochemistry analysis of paraffi...



Western blot analysis of lysates from He...