

Product Datasheet

PAK4 (phospho-Ser474) antibody (orb767332)

Description Rabbit polyclonal antibody to PAK4 (phospho-Ser474)

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, IHC-P, WB

Immunogen The antiserum was produced against synthesized peptide derived from human PAK4/5/6 around the phosphorylation site of Ser474. AA range:441-490

Preservatives PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Concentration 1 mg/ml

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.


Note For research use only

Isotype IgG


Clonality Polyclonal

MW 64

Expiration Date 12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Western Blot analysis of various, using ...