



Product Datasheet

CPTI-M antibody (orb767484)


Description	Rabbit polyclonal antibody to CPTI-M
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human CPT1B. AA range:401-450
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	88
Expiration Date	12 months from date of receipt.



Western Blot analysis of HepG2 cells using...



Western Blot analysis of NIH-3T3 cells using...



Western blot analysis of CPT1B Antibody....