

Product Datasheet

Pseudomonas aeruginosa serotype 5c antibody (orb76792)

Description	Mouse monoclonal antibody to Pseudomonas aeruginosa serotype 5c
Species/Host	Mouse
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, ICC
Immunogen	Native Pseudomonas aeruginosa serotype 5c
Target	Pseudomonas aeruginosa serotype 5c
Preservatives	Purified antibody containing PBS 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	SP2/0
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	Optimal antibody dilution should be determined by user laboratory since it may be dependent on the detection system used in both ELISA and IF
Expiration Date	12 months from date of receipt.