

## Product Datasheet

# Pseudomonas aeruginosa serotype 6c Monoclonal Antibody [MJ6302] (orb76793)

<b>Description</b>	Mouse monoclonal antibody to Pseudomonas aeruginosa serotype 6c
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Bacteria
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IF
<b>Immunogen</b>	Native Pseudomonas aeruginosa serotype 6c
<b>Target</b>	Pseudomonas aeruginosa serotype 6c
<b>Preservatives</b>	Purified antibody containing PBS 0.1% sodium azide
<b>Concentration</b>	1mg/ml
<b>Note</b>	For research use only
<b>Application notes</b>	SP2/0
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Clone Number</b>	MJ6302
<b>Dilution Range</b>	Optimal antibody dilution should be determined by user laboratory since it may be dependent on the detection system used in both ELISA and IF
<b>Expiration Date</b>	12 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)