


Product Datasheet

SUN2 antibody (orb76819)

Description	Rabbit monoclonal antibody to SUN2
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Mouse SUN2 residues 1-18
Target	SUN2
Preservatives	Serum containing 0.1% sodium azide
Note	For research use only
Clonality	Monoclonal
Uniprot ID	Q9UH99
Entrez	25777
Dilution Range	Optimal antibody dilution should be determined by titration, however, as a guideline try 1:1000 for immunoblotting and 1:500 for immunofluorescence
Expiration Date	12 months from date of receipt.



Confocal immunofluorescence analysis of ...