



Product Datasheet

IL15R alpha antibody (orb768774)


Description	Rabbit polyclonal antibody to IL15R alpha
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human IL15RA. AA range:99-148
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	32
Expiration Date	12 months from date of receipt.



Western Blot analysis of COLO cells using...



Western blot analysis of various cell lysates using...



Western blot analysis of lysates from COLO cells using...