



Product Datasheet

NT5C3 antibody (orb769381)


| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to NT5C3 |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IF, IHC-P, WB |
| Immunogen | The antiserum was produced against synthesized peptide derived from human NT5C3. AA range:11-60 |
| Preservatives | PBS with 0.02% sodium azide and 50% glycerol pH 7.4. |
| Concentration | 1 mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 38 |
| Expiration Date | 12 months from date of receipt. |



Immunohistochemistry analysis of paraffi...



Western blot analysis of lysates from Ju...



Western blot analysis of the lysates fro...