



Product Datasheet

PDGFR-beta (phospho-Tyr771) antibody (orb769416)

Description	Rabbit polyclonal antibody to PDGFR-beta (phospho-Tyr771)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized phospho-peptide around the phosphorylation site of human PDGFR- β (phospho Tyr771)
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	125
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Western Blot analysis of mouse lung, rat...