



Product Datasheet

PI 3-kinase p85 beta (phospho-Tyr464) antibody (orb769467)


Description	Rabbit polyclonal antibody to PI 3-kinase p85 beta (phospho-Tyr464)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized phospho-peptide around the phosphorylation site of human PI 3-kinase p85β (phospho Tyr464)
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	85
Expiration Date	12 months from date of receipt.



Enzyme-Linked Immunosorbent Assay (Phosp...



Western blot analysis of lysate from Jur...



Immunohistochemical analysis of paraffin...