



Product Datasheet

Axl (phospho-Tyr691) antibody (orb769645)


Description	Rabbit polyclonal antibody to Axl (phospho-Tyr691)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human AXL around the phosphorylation site of Tyr691. AA range:657-706
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	110 97
Expiration Date	12 months from date of receipt.



Western Blot analysis of HuvEc cells treated with HUVEC EGF 200ng/ml 15'



Immunohistochemistry analysis of paraffin-embedded tissue sections



Western blot analysis of lysates from HUVEC cells