



Product Datasheet

PKC delta (phospho-Tyr64) antibody (orb769648)


Description	Rabbit polyclonal antibody to PKC delta (phospho-Tyr64)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human PKCD around the phosphorylation site of Tyr64. AA range:30-79
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	77
Expiration Date	12 months from date of receipt.



Enzyme-Linked
Immunosorbent
Assay (Phospho...



Western blot analysis
of lysates from He...



Immunohistochemical
analysis of paraffin...