

Product Datasheet

PKD1 (phospho-Tyr463) antibody (orb769662)

Description Rabbit polyclonal antibody to PKD1 (phospho-Tyr463)

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, IHC-P, WB

Immunogen The antiserum was produced against synthesized peptide derived from human PKD1/PKC mu around the phosphorylation site of Tyr463. AA range:429-478

Preservatives PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Concentration 1 mg/ml

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.


Note For research use only

Isotype IgG


Clonality Polyclonal

MW 130


Expiration Date 12 months from date of receipt.



Enzyme-Linked
Immunosorbent Assay
(Phosp...



Immunohistochemistry
analysis of paraffi...



Western blot analysis
of lysates from He...