



Product Datasheet

DNA-PKCS (phospho-Thr2647) antibody (orb769678)


Description	Rabbit polyclonal antibody to DNA-PKCS (phospho-Thr2647)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P
Immunogen	The antiserum was produced against synthesized peptide derived from human DNA-PK around the phosphorylation site of Thr2647. AA range:2613-2662
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.



Enzyme-Linked
Immunosorbent Assay
(Phosp...



Immunofluorescence
analysis of HUVEC
cel...



Immunohistochemistry
analysis of paraffi...