



Product Datasheet

BAF250b antibody (orb769784)


Description	Rabbit polyclonal antibody to BAF250b
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human BAF250B. AA range:1371-1420
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	170
Expiration Date	12 months from date of receipt.



Western Blot analysis of 3T3 cells using...



Immunohistochemical analysis of paraffin...



Immunofluorescence analysis of HUVEC cel...