



Product Datasheet

PSK2 antibody (orb769794)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to PSK2 |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IHC-P, WB |
| Immunogen | The antiserum was produced against synthesized peptide derived from human TAOK1. AA range:431-480 |
| Preservatives | PBS with 0.02% sodium azide and 50% glycerol pH 7.4. |
| Concentration | 1 mg/ml |
| Storage | Store at -20°C. Avoid repeated freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 120 |
| Expiration Date | 12 months from date of receipt. |



Immunohistochemistry analysis of paraffi...



Western blot analysis of lysates from NI...