



Product Datasheet

GABA T-1 antibody (orb770093)

Description	Rabbit polyclonal antibody to GABA T-1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human SLC6A1. AA range:381-430
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	67
Expiration Date	12 months from date of receipt.



Western blot analysis of SLC6A1 Antibody...



Western blot analysis of the lysates fro...