


Product Datasheet

ADRA1D antibody (orb77040)

Description	Rabbit polyclonal antibody to ADRA1D
Species/Host	Rabbit
Reactivity	Drosophila, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	A peptide corresponding to middle region of human ADRA1D protein
Target	ADRA1D
Form/Appearance	Lyophilized powder with 0.02% NaN ₃ .
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Reconstitute in 100µL of sterile H ₂ O with 50% glycerol.
Isotype	IgG
Clonality	Polyclonal
MW	60kDa
Dilution Range	WB: 1:200-1:1,000
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 and MCF-7 ...